# Nettundervisning

## Dekontaminering - fleksible endoskop

Linda Ashurst  
Seksjon for dekontaminering  
Avd. for smittevern

# Disposisjon

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Faser i dekontamineringsprosessen

Testing av vaskedekontaminator iht. NS-EN 15883-4

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
# Vannbehandlingsystemer for sluttskyllefasen (UK, HTM 01-06)

Figure 6 Selection of water treatment systems for EWD final rinse

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
# Proessen i vaskedekontaminator

# NS-EN ISO 15883: 2009

**5.2.1.1** During at least part of each of the cleaning, disinfection and rinsing phases, the device channel irrigation system shall ensure that the various process fluids flow through each of the internal channels and/or cavities of the devices that are required to be cleaned and disinfected. Assurance that this has taken place shall be provided either by:

 a) the automatic controller providing means to verify the flow of process fluids through each channel (see 5.2.2);

or

 b) requiring in the instructions for use that the user:

- verifies that all channels allow the free passage of water before the device is loaded into the WD;
- confirms that all necessary connections were made before, and were still in place at the end of, the cycle;
- confirms by reference to the WD process record that the supply of process fluids was maintained during each stage of the process (see 5.5):
- verifies flow through each endoscope channel at the end of each operational cycle or immediately before use.

 Option a) above is the preferred method. Before choosing equipment conforming with option b) users should consider the increased requirement for staff training, staff time during reprocessing and the lack of an independent record that the process was carried out in a satisfactory manner.

# NS-EN ISO 15883: 2009

## 5.2.2 Verification of device channel irrigation by the automatic controller

**5.2.2.1** The WD manufacturer shall specify for each channel the maximum extent of flow reduction permissible (e.g. change in flow volumes, pressures, rates, etc.) that will not impair the efficacy of the process [see 8 a)].

The WD manufacturer should request the device manufacturer to supply relevant data e.g. dimensions of connectors, internal dimensions of channels, maximum pressures to which channels may be subjected, in order to enable the WD manufacturer to determine the flow that will occur through the unobstructed channels.

When one or more channels of the device are obstructed to an extent that would impair the efficacy of the process, the automatic controller shall cause a fault to be indicated.

Compliance with this requirement shall be demonstrated by testing in accordance with 6.6.

With some designs of endoscope a blockage in one channel may cause the flow to be diverted to another channel or port. Under these circumstances detection of an obstruction by the automatic controller may not be reliable. The user should refer to the device manufacturer's instructions for the method to be used to verify that all channels are free from obstructions.

**5.2.2.2** When one or more channels of the device are not connected to the WD, the automatic controller shall cause a fault to be indicated.

Compliance with this requirement shall be demonstrated by testing in accordance with 6.7.

**5.2.2.3** The automatic controller shall verify that the duration of flow of the relevant process fluids met or exceeded the minimum exposure times established during validation studies as necessary for each process stage. Failure to achieve the required flow shall cause a fault to be indicated.

**5.2.2.4** When there is a common connection for fluid at the same supply pressure to more than one channel, the WD manufacturer shall provide evidence that the flow through each of the channels meets, or exceeds, the minimum required for effective cleaning, disinfection and rinsing of each device channel to be processed.

# NS-EN ISO 15883: 2009

## 5.5 Process chemicals

The conditions of use (temperature, concentration etc.) within the WD for all process chemicals (detergent, disinfectant etc.) shall be within the limits specified by the process chemical manufacturer.

Where required for testing purposes the WD manufacturer shall obtain from the process chemical manufacturer the method(s) to be used to neutralize the process chemical (e.g. to stop any further anti-microbial activity) and shall make this information available to the user (see 4.4.2.3, NOTE).

## 5.6 Process verification

Process verification shall be in accordance with ISO 15883-1:2006, 5.11.4 c).

## 5.7 Dosing systems

For those WDs in which the required dose of process chemical is contained in a single-dose container which is replaced before each cycle, means shall be provided to ensure that the intended volume has been dispensed.

Compliance shall be tested in accordance with 6.10.

For WDs in which process chemicals are supplied in multi-dose containers, ISO 15883-1:2006, 5.7 applies.

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