

Desinfeksjonsmidler som anbefales til ultralydprober og som produsentene har bekreftet kan benyttes uten å skade deres prober. November 2024			Desinfeksjonsmidler til teknisk bruk <i>Bruksområder, effekt og virketider for desinfeksjonsmidler</i> Infeksjonskontroll.no og <i>desinfeksjonskluter</i> Infeksjonskontroll.no		
Produsent av ultralydprober	Varemerker	Leverandør	Pereddiksyre	Klordioksid	Etanol brukes bare under rene forhold
Siemens Healthineers	ACUSON Sequoia Ultrasound System	Siemens Healthineers	Rely+On PeraSafe	Tristel Trio Wipe, Tristel Duo ULT	
BK Medical	bkSpecto	KeboMed AS	Rely+On PeraSafe	Tristel Trio Wipe, Tristel Duo ULT	Desinfeksjonssprit 70 - 75 % <i>kun til ikke kritisk utstyr</i>
FUJIFILM SonoSite	SonoSite	Secma AS		Tristel Trio Wipe, Tristel Duo ULT	
GE Healthcare, GE Vingmed Ultrasound	GE Healthcare	GE Healthcare Norge AS, GE Vingmed Ultrasound	Rely+On PeraSafe	Tristel Trio Wipe, Tristel Duo ULT	Desinfeksjonssprit 70 - 75 % <i>kun til ikke kritisk utstyr</i>
Philips	Philips ATL	Philips Norge AS		Tristel Trio Wipe, Tristel Duo ULT	
Hitachi – Aloka	Aloka	Scan-Med AS Norway	Rely+On PeraSafe	Tristel Trio Wipe, Tristel Duo ULT	
Canon Medical	Canon Aplio, Xario, Aviamo	Tromp Medical AS	Rely+On PeraSafe, Clinell Peracetic Acid Wipes	Tristel Trio Wipe, Tristel Duo ULT	Desinfeksjonssprit 70 - 75 % <i>kun til ikke kritisk utstyr</i>
Samsung Healthcare	Samsung	Inter Medical AS		Tristel Trio Wipe, Tristel Duo ULT LifeClean	
Natus Neuro	Natus® InVisus™ Ultrasound Sonara	Timik AS		Tristel Trio Wipe, Tristel Duo ULT	